

Micropump MPmlh+

Reference	KP1.12.220.6 MP+ mlh Multi syringes KP1.12.202.4 MP+ mlh BD 20 ml KP1.12.209.4 MP+ mlh BBraun 20 ml KP1.12.208.4 MP+ mlh Codan 20 ml
Delivery modes	- Continuous - Continuous + bolus
Drive Accuracy	+2%
Safety functions	- Multi point(s) syringe dislodged detection - Lockable flow rate
Downloadable event log	Yes
Specifications	For MPmlh+ Multi syringes: uses manufacturer's standard syringes (BD, Terumo, BBraun, Codan/LMT) from 10 to 50 mL For the rest MPmlh+ variants: each variant dedicated to specific syringes • Fail-safe dual micro-processor motor monitoring • Backlit LCD confirms correct operation with alarm and warning indicators
Alarms & warnings	- Pre-end of infusion - Low battery - Occlusion - System failure - End of infusion - Depleted battery - Pump on hold - Pump tamper warning - Priming indication - Bolus delivery
Rate	0.1 – 28.0 ml/hr, in 0.1 ml/hr steps (for BD-20 ml)
Battery life	More than 70 full infusions at intermediate rate (2 to 3 months continuous)
Weight with battery	220 g including batteries (Alkaline LR03 AAx6)
Dimensions	170x61x32 mm

Includes 6 x 1,5V alkaline batteries & carrying bag

Accessories for MPmlh+ Pump

	<p>MP+ Carrying bag with transparent window 20-50 ml Included in all MP+ Micropump systems</p> <p>KM1.YY.604.1</p>		<p>Disposable bag for MP+ pumps (Box of 20 pcs)</p> <p>KM1.YY.561.1</p>	 <p>MP+ Soft Locker box 20 ml</p> <p>KS1.12.164.1</p>
	<p>MP+ Pole Clamp</p> <p>KM1.EE.121.1</p>			

96083 Rev.3-Rel.Date: 2019.10.29